

Medcare Transfer & Position Aids

MEDCARE
PRODUCTS

Medsheet & Medsheet Tube

The Medsheet is composed of two nylon sheets, which allow for easy movement of a patient during repositioning.

Item Number:

400068: Medsheet set of two
Dimensions: 34" W x 67" L

The Medsheet Tube is a loop of fabric with webbing handles used for repositioning a patient in a bed.

Item Number:

400076: Medsheet Tube
Dimensions: 32" W x 40" L
400095: Medsheet Tube
Dimensions: 34" W x 67" L




One-Way-Glide Anti-Slip Sheets

The One-Way-Glide is used to help patients maintain good position in the wheelchair/day chair/bed. Available in 3 sizes, it has been designed to move in the direction of the arrow only, and helps prevent the patient from slipping forward.

Item Number:

255200: Small chair size
Dimensions: 14.5" W x 17" L
255201: Medium chair size
Dimensions: 17.75" W x 17.75" L
255202: Large chair size
Dimensions: 23.5" W x 31.5" L



Gait/Handling Belt

A Gait or Handling Belt is an ideal aid when assisting with transfers using transfer boards or turning discs and can be used with one or two caregivers. It enables better handling techniques while assisting in transfers or repositioning.

Item Number:

402555: Medium Gait Belt
Waist Size: 34"-49"
402560: Large Gait Belt
Waist Size: 42"-60"



Gait Belt Features:

- Non-skid inner material that has superior hold
- Angled handles to ensure quick grab of transferee
- Ripstop technology on outer shell
- Extra wide belt provides greater contact area with patient
- Inner comfort and support padding

“ At Medcare Products, our mission is to provide the highest-quality patient lifting, transferring, and positioning equipment, along with the highest standard of training and service in the industry.”



Transfer & Position Aids (cont.)

Motility Board

The Motility Board is ideal for assisted and non-assisted transfers and allows easy maneuvering from bed to chair, chair to wheelchair, etc. It also has anti-slip pads and has a built in antimicrobial agent for increased infection control. The maximum safe working load is 600 lbs.

Item Number:

255701: Motility Board




Lateral Assist

The Lateral Assist is designed for supine lateral transfers from bed to bed/stretcher/theatre table/x-ray/shower trolley. The collapsible board is made of plasticized inner board that is sealed with electronic stitching, making it water resistant. Each board comes with a nylon cover. Cleaning can be done with a 0.5% accelerated hydrogen peroxide.

Item Number:

100300: 2 Move Boards with Nylon Cover
Dimensions: 15.5" W x 63" L



Disposable Slide Sheets

Disposable Slide Sheets are designed for multi-use with the same patient and then disposed of when soiled or no longer required.

Item Number:

600310: Pack of 100
Dimensions: 39" W x 78" L



Multiglide Glove

The Multiglide Glove is an effective aid in reducing friction when handling patients. The glove is carefully glided under the patient and is used for repositioning or for checking therapeutic seating surfaces.

Item Number:

600309: Pair of Gloves




Onward Turntable

The Onward Turntable is a non-slip, injection molded turning disc. It is positioned on the floor as an aid to turn a weight bearing person. The disc is easy to disinfect between uses.

Item Number:

250006: Onward Turntable
Diameter: 13.25"



Medcare Products

151 East Cliff Road
Burnsville, MN 55337

Toll-free: 800 695 4479
Fax: 952 300 2100

www.MedcareLifts.com

© 2014 Medcare Products Inc.

Rev 031914

IMPORTANT: Always assess the patient to ensure the correct aid is being used. Transfer and reposition aids are always at the discretion of the facility. For tips and training on using Medcare's various transfer and position aids, contact your local Safe Patient Handling Consultant or Medcare's technical support at 1-800-695-4479.